

Product information

Compression Bandage

The absolutely uncomplicated, ready-to-use SÖHNGEN® pressure bandage for all injuries with flowing bleeding. The special pressure bandage pack replaces any conventional pressure bandage, it does not slip or shift and can be applied quickly and correctly even by inexperienced users. The special pressure bandage pack has a pressure pad with a firm, insulated core and very strong absorbent layers, which prevents bleeding through. The sewn-on elastic WS fixation bandage is extremely secure and non-slip without impeding blood circulation. The special pressure bandage pack ensures that bleeding stops perfectly, there is no congestion and no strangulation. A safe, reliable and successful method without time-consuming preparation. Colour code orange, EO-sterilized, individually packed in peel packs.



Art.-No	1002001
GTIN	4250108803157
Custom Tariff No	30059050
Country of origin	DE-HE
eclass 4	34030301
eclass 5	34030101
Height in mm	160
Width in mm	110
Depth in mm	110
Weight in kg	0,042
Norm	
Colour	
Material	
Water resistant	
Lockable	no
Splashproof	
Wall mount	none

Drug	no drug
Medical device class	MD class I(s)
Biocide	no
PPE	
PZN	03147548
BfArM appr. No	
BAuA Reg. No	
EMDN	M0399
Base UDI	42501088009J4
UDI-DI	04250108803157
UN-No	
O2 in liters	0
Manufacturer	W.Söhnngen GmbH
Manufacturer Website	www.soehngen.com
Manufacturer Contact	info@soehngen.com