

Product information

Dressing Sheet SO 60 x 80 cm

Dressing cloth SO. PP non-woven, absorbency approx. 125 g/m². Meets the minimum requirements of DIN 13 152. Individually sterile in peel pack.



Art.-No	1003009
GTIN	4250108808138
Custom Tariff No	30059050
Country of origin	DE-HE
eclass 4	34030301
eclass 5	34030101
Height in mm	159
Width in mm	125
Depth in mm	11
Weight in kg	0,018
Norm	DIN 13152
Colour	
Material	
Water resistant	
Lockable	no
Splashproof	
Wall mount	none

Drug	no drug
Medical device class	MD class I(s)
Biocide	no
PPE	
PZN	
BfArM appr. No	
BAuA Reg. No	
EMDN	M040299
Base UDI	42501088006HW
UDI-DI	04250108808138
UN-No	
O2 in liters	0
Manufacturer	W.Söhngen GmbH
Manufacturer Website	www.soehngen.com
Manufacturer Contact	info@soehngen.com