

Product information

aluderm® Dressing Sheet 80 x 120cm

Aluminum vaporized, individually sterile in peel pack. aluderm® - the gentle wound therapy. aluderm® with aluminum-vaporized wound dressing is ideal for all superficial wounds and injuries, superficial burns, abrasions, lacerations and cuts. aluderm® is not just a brand name, it is the epitome of perfect wound care. Whether minor injuries or larger surface wounds - aluderm® always brings relief. Medical professionals and patients are enthusiastic about aluderm®. The success of aluderm® in use in specialist clinics confirms the exceptional nature and quality of this dressing material. aluderm® is an extraordinarily soft and supple dressing material: aluminium-vaporized wound pad made of polyester nonwoven, very strong absorbent layer (absorbency 1,400 g/m²) and hydrophobic PP final layer. aluderm® has many good properties: aluderm® does not stick to the wound, does not fray, is highly absorbent, breathable, hypoallergenic, can help with skin formation and healing, physiologically harmless, without optical brighteners and free of chemical binders.



Art.-No	1003206
GTIN	4250108803225
Custom Tariff No	30059050
Country of origin	DE-HE
eclass 4	34030301
eclass 5	34030101
Height in mm	175
Width in mm	120
Depth in mm	50
Weight in kg	0,132
Norm	DIN 13152
Colour	white
Material	
Water resistant	
Lockable	no
Splashproof	
Wall mount	none

Drug	no drug
Medical device class	MD class I(s)
Biocide	no
PPE	
PZN	03147353
BfArM appr. No	
BAuA Reg. No	
EMDN	M040202
Base UDI	42501088007HY
UDI-DI	04250108803225
UN-No	
O2 in liters	0
Manufacturer	W.Söhngen GmbH
Manufacturer Website	www.soehngen.com
Manufacturer Contact	info@soehngen.com