

Product information

aluderm®aluplast Wound Plaster 5 m x 8 cm, elastic

aluderm®-aluplast wound dressing plaster - rolled. 5 m plaster strips in one piece, rolled up, unsealed. aluderm®-aluplast is the special wound dressing plaster with aluminum-vaporized aluderm® wound pad and absorbent fleece. aluderm® does not stick to the wound, does not fray, is highly absorbent, can help with skin formation and healing, is breathable. aluderm®-aluplast elastic: Backing material highly elastic polyester fleece, white, particularly skin-friendly, good air permeability, water-repellent. The polyacrylate adhesive is additionally activated by body heat. Particularly supple, adapts to all movements, hypoallergenic, non-irritant. Also particularly suitable for people with sensitive skin or plaster allergies.



Art.-No	1009190
GTIN	4250108808794
Custom Tariff No	30051000
Country of origin	DE-HE
eclass 4	34030301
eclass 5	34030101
Height in mm	5000
Width in mm	3
Depth in mm	80
Weight in kg	0,114
Norm	DIN 13019
Colour	white
Material	
Water resistant	
Lockable	no
Splashproof	
Wall mount	none

Drug	no drug
Medical device class	MD class I
Biocide	no
PPE	
PZN	00007841
BfArM appr. No	
BAuA Reg. No	
EMDN	M04010101
Base UDI	42501088113HY
UDI-DI	04250108808794
UN-No	
O2 in liters	0
Manufacturer	W.Söhngen GmbH
Manufacturer Website	www.soehngen.com
Manufacturer Contact	info@soehngen.com